

Light Adapter Guide

Revision Date: 04/17/2015










Two types of lights:

- Dove Tail LED lights for TTL loupes
- Clip LED/Fiber Optic lights for FLM loupes

Standard adapters for TTL loupes				
Part no.	Photo		Description call-out	Used on
25794			REVISED Aero Dove Tail Adapter (1)	SurgiTel Aero SurgiTel ComfortFlex Oakley Deringer
N/A			NEW Aero Dove Tail Adapter (2)	SurgiTel Aero (included with TTL purchase, frame modification required)
25798-1	 <i>Screw type 1</i>		Ergo Dove Tail Adapter	SurgiTel Ergo
25798-2	 <i>Screw type 2</i>		Orascoptic Rudy Frames Dove Tail Adapter	Orascoptic Rudy Frames
25798-3	 <i>Screw type 3</i>		Orascoptic Metal Frames Dove Tail Adapter	Orascoptic Metal Frames
25798-4	 <i>Screw type 4</i>		Designs for Vision Wire Frames Dove Tail Adapter	DFV Wire Frames
25798-5	 <i>Screw type 5 (modified adapter)</i>		NEW Zeiss FLM Loupe Dove Tail Adapter	Zeiss FLM Loupe
25793			REVISED Half Jacket/Flak Jacket/ Prime Dove Tail Frame Adapter	Oakley Flak Jacket Oakley Half Jacket (2.0) Sundog Prime

Standard adapters for TTL loupes (continued)				
Part no.	Photo		Description call-out	Used on
25795			Slide-retainer Dove Tail Adapter	Oakley Radar
25864			Exam Vision Dove Tail Adapter	Exam Vision Frames (all)
25865			Merident FLM Loupes Dove Tail Adapter	Merident FLM Loupes
25932			Designs for Vision Nike Skylon Dove Tail Adapter	DFV Nike Skylon Frame
25169			SurgiCam Adapter to Fixed Headband Light with Iris	SurgiCam
Clip adapters for SurgiTel TTL loupes (available per request)				
25492			Aero Clip Adapter	SurgiTel Aero SurgiTel ComfortFlex Oakley Deringer
25207			Half Jacket/Flak Jacket Clip Adapter	Oakley Flak Jacket Oakley Half Jacket (2.0)
23815-1			Ergo Slide-retainer Clip Adapter	SurgiTel Ergo

Clip adapters for SurgiTel TTL loupes (available per request) (continued)				
Part no.	Photo		Description call-out	Used on
23815-2			Radar Slide-retainer Clip Adapter	Oakley Radar All frames with slide-retainer
<i>Note for part no. 23815-X: These two adapters are not the same. Part no. 23815-1 includes an Ergo slide-retainer and part no. 23815-2 does not.</i>				
Competitor Adapters				
24741			Designs for Vision Yeoman Frame Dove Tail & Clip Adapter	DFV Yeoman Frame
24742			Orascoptic Titanium Frames Clip Adapter	Orascoptic Titanium Frames
24739			Orascoptic FLM Loupes Clip Adapter	Orascoptic FLM Loupes
25216			Heine FLM Loupes Clip Adapter	Heine FLM Loupes
24740			Keeler FLM Loupes Clip Adapter	Keeler FLM Loupes
25217			Zeiss FLM Loupes Clip Adapter	Zeiss FLM Loupes

Universal Light Adapters				
Part no.	Photo		Description call-out	Used on
23800-2			TTL Clip Adapter	DFV Wire Frames Merident Frames (all) Zeiss Frames (all) Exam Vision Frames (all) All frame types (majority)
23800-4			TTL Sturdy Dove Tail Adapter	DFV Wire Frames Merident Frames (all) Zeiss Frames (all) Exam Vision Frames (all) All frame types (majority)
25805			TTL Lightweight Dove Tail Adapter	DFV Yeoman Frame Merident Frames (all) Zeiss Frames (all) All frame types (majority)
Safety Frame				
Part no.	Photo		Description call-out	Used on
25874			Smith & Wesson Frame with Dove Tail Adapter	Smith & Wesson Frame
Replacement LED Light Unit Adapters				
25975			Dove Tail Light Unit Adapter	SurgiTel LED Light Unit
26034			Clip Light Unit Adapter	SurgiTel LED Light Unit